

Antibodies Product Datasheet

Product Name:	p97 MAPK Antibody
Catalog Number:	R12-3293
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	p97 MAPK Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human p97 MAPK.
Concentration:	1 mg/ml
Gene Name:	MAPK6
Other Name:	ERK-3; ERK3; Extracellular signal- regulated kinase 3; MAP kinase isoform p97; MAPK6; MK06; Mitogen-activated protein kinase 6; PRKM4; PRKM6; kinase ERK3; p55-MAPK; p97-MAPK
Human Swiss Prot No.:	Q16659
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	p97 MAPK Antibody antibody detects endogenous levels of p97 MAPK.
Storage and Stability:	Stable at -20°C for at least 1 year.