

## Antibodies Product Datasheet

<b>Product Name:</b>	Brg-1 Monoclonal Antibody
<b>Catalog Number:</b>	N0670
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant human Brg-1 protein fragments expressed in E.coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	SMARCA4
<b>Protein Name:</b>	Transcription activator BRG1
<b>Other Name:</b>	SMARCA4; BAF190A; BRG1; SNF2B; SNF2L4; Transcription activator BRG1; ATP-dependent helicase SMARCA4; BRG1-associated factor 190A; BAF190A; Mitotic growth and transcription activator; Protein BRG-1; Protein brahma homolog 1; SNF2-beta; SWI/S
<b>Human Swiss Prot No.:</b>	P51532
<b>Formulation:</b>	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
<b>Specificity:</b>	Brg-1 Monoclonal Antibody detects endogenous levels of Brg-1 protein.
<b>Storage and Stability:</b>	-20°C/1 year