
Antibodies Product Datasheet

Product Name:	EP1 Antibody
Catalog Number:	C30225
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human PTGER1. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	PTGER1
Protein Name:	Prostaglandin E2 receptor EP1 subtype
Other Name:	PTGER1; Prostaglandin E2 receptor EP1 subtype; PGE receptor EP1 subtype; PGE2 receptor EP1 subtype; Prostanoid EP1 receptor
Human Swiss Prot No.:	P34995
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EP1 Polyclonal Antibody detects endogenous levels of EP1 protein.
Storage and Stability:	-20°C/1 year