

Antibodies Product Datasheet

Product Name:	LOX12 Antibody
Catalog Number:	C51161
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal LOX12
Concentration:	1 mg/ml
Gene Name:	ALOX12 LOG12
Protein Name:	Arachidonate 12-lipoxygenase, 12S-type (12S-LOX) (12S-lipoxygenase) (EC 1.13.11.31) (Platelet-type lipoxygenase 12)
Other Name:	Arachidonate 12-lipoxygenase, 12S-type (12S-LOX) (12S-lipoxygenase) (EC 1.13.11.31) (Platelet-type lipoxygenase 12)
Human Swiss Prot No.:	P18054
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LOX12 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year