

Antibodies Product Datasheet

Product Name:	FoxO1/3 Antibody
Catalog Number:	B0485
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FOXO1A/3A. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	FOXO1/FOXO3
Protein Name:	Forkhead box protein O1/3
Other Name:	FOXO1; FKHR; FOXO1A; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma; FOXO3; FKHL1; FOXO3A; Forkhead box protein O3; AF6q21 protein; Forkhead in rhabdomyosarcoma-like 1
Human Swiss Prot No.:	Q12778/O43524
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FoxO1/3 Polyclonal Antibody detects endogenous levels of FoxO1/3 protein.
Storage and Stability:	-20°C/1 year