

Antibodies Product Datasheet

Product Name:	KCNF1 Antibody
Catalog Number:	C17808
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KCNF1. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	KCNF1
Protein Name:	Potassium voltage-gated channel subfamily F member 1
Other Name:	KCNF1; Potassium voltage-gated channel subfamily F member 1; Voltage-gated potassium channel subunit Kv5.1; kH1
Human Swiss Prot No.:	Q9H3M0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KCNF1 Polyclonal Antibody detects endogenous levels of KCNF1 protein.
Storage and Stability:	-20°C/1 year