

Antibodies Product Datasheet

Product Name:	BART1 Antibody
Catalog Number:	C14350
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ARL2BP. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	ARL2BP
Protein Name:	ADP-ribosylation factor-like protein 2-binding protein
Other Name:	ARL2BP; BART; BART1; ADP-ribosylation factor-like protein 2-binding protein; ARF-like 2-binding protein; Binder of ARF2 protein 1
Human Swiss Prot No.:	Q9Y2Y0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BART1 Polyclonal Antibody detects endogenous levels of BART1 protein.
Storage and Stability:	-20°C/1 year