

Antibodies Product Datasheet

Product Name:	Ub (Acetyl-Lys63) Antibody
Catalog Number:	D12236
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Ub (Acetyl Lys63)
Concentration:	1 mg/ml
Gene Name:	UBA52 UBCEP2
Protein Name:	Ub (Acetyl Lys63)
Other Name:	Ubiquitin-60S ribosomal protein L40 (CEP52) (Ubiquitin A-52 residue ribosomal protein fusion product 1) [Cleaved into: Ubiquitin; 60S ribosomal protein L40]
Human Swiss Prot No.:	P62987/P62979/P0CG47/P0CG48
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat Ub (Acetyl Lys63)
Storage and Stability:	-20°C/1 year