

Antibodies Product Datasheet

Product Name:	SAPS3 Antibody
Catalog Number:	C49190
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SAPS3
Concentration:	1 mg/ml
Gene Name:	PPP6R3 C11orf23 KIAA1558 PP6R3 SAPL SAPS3
Protein Name:	Serine/threonine-protein phosphatase 6 regulatory subunit 3 (SAPS domain family member 3) (Sporulation-induced transcript 4-associated protein SAPL)
Other Name:	
Human Swiss Prot No.:	Q5H9R7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SAPS3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year