

## Antibodies Product Datasheet

<b>Product Name:</b>	PKM2 Antibody
<b>Catalog Number:</b>	C30826
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC-p, IF(paraffin section), WB
<b>Dilutions:</b>	IHC-p 1:50-200, WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human PKM2
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PKM OIP3 PK2 PK3 PKM2
<b>Protein Name:</b>	PKM2
<b>Other Name:</b>	Pyruvate kinase isozymes M1/M2 (EC 2.7.1.40) (Cytosolic thyroid hormone-binding protein) (CTHBP) (Opa-interacting protein 3) (OIP-3) (Pyruvate kinase 2/3) (Pyruvate kinase muscle isozyme) (Thyroid hormone-binding protein 1) (THBP1) (Tumor M2-PK) (p58)
<b>Human Swiss Prot No.:</b>	P14618
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of human PKM2
<b>Storage and Stability:</b>	-20°C/1 year