
Antibodies Product Datasheet

Product Name:	Dmc1 Antibody
Catalog Number:	C13043
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DMC1. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	DMC1
Protein Name:	Meiotic recombination protein DMC1/LIM15 homolog
Other Name:	DMC1; DMC1H; LIM15; Meiotic recombination protein DMC1/LIM15 homolog
Human Swiss Prot No.:	Q14565
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Dmc1 Polyclonal Antibody detects endogenous levels of Dmc1 protein.
Storage and Stability:	-20°C/1 year