

Antibodies Product Datasheet

Product Name:	MAGE1 Antibody
Catalog Number:	C44659
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 760-840
Concentration:	1 mg/ml
Gene Name:	MAGEE1 HCA1 KIAA1587
Protein Name:	Melanoma-associated antigen E1 (Alpha-dystrobrevin-associated MAGE Protein) (DAMAGE) (Hepatocellular carcinoma-associated protein 1) (MAGE-E1 antigen)
Other Name:	
Human Swiss Prot No.:	Q9HCl5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	MAGE1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year