

Antibodies Product Datasheet

Product Name:	Gab 1 Antibody
Catalog Number:	B0788
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GAB1. AA range:644-674
Concentration:	1 mg/ml
Gene Name:	GAB1
Protein Name:	GRB2-associated-binding protein 1
Other Name:	GAB1; GRB2-associated-binding protein 1; GRB2-associated binder 1; Growth factor receptor bound protein 2-associated protein 1
Human Swiss Prot No.:	Q13480
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Gab 1 Polyclonal Antibody detects endogenous levels of Gab 1 protein.
Storage and Stability:	-20°C/1 year