

Antibodies Product Datasheet

Product Name:	MRPL51 Antibody
Catalog Number:	R12-3045
Reactivity:	Human, Mouse
Applications:	WB, IHC, ELISA
Dilutions:	MRPL51 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRPL51.
Concentration:	1 mg/ml
Gene Name:	MRPL51
Other Name:	39S ribosomal protein L51 mitochondrial; L51mt; MRP-L51; bMRP64; bMRP-64
Human Swiss Prot No.:	Q4U2R6
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	MRPL51 Antibodyantibody detects endogenous levels of MRPL51.
Storage and Stability:	Stable at -20°C for at least 1 year.