

Antibodies Product Datasheet

Product Name:	PPAPDC3 Antibody
Catalog Number:	C48983
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PPAPDC3
Concentration:	1 mg/ml
Gene Name:	PPAPDC3 C9orf67
Protein Name:	Probable lipid phosphate phosphatase PPAPDC3 (EC 3.1.3.-) (Phosphatidic acid phosphatase type 2 domain-containing protein 3)
Other Name:	
Human Swiss Prot No.:	Q8NBV4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PPAPDC3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year