

Antibodies Product Datasheet

Product Name:	Sin3B Antibody
Catalog Number:	C10163
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SIN3B. AA range:221-270
Concentration:	1 mg/ml
Gene Name:	SIN3B
Protein Name:	Paired amphipathic helix protein Sin3b
Other Name:	SIN3B; KIAA0700; Paired amphipathic helix protein Sin3b; Histone deacetylase complex subunit Sin3b; Transcriptional corepressor Sin3b
Human Swiss Prot No.:	O75182
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Sin3B Polyclonal Antibody detects endogenous levels of Sin3B protein.
Storage and Stability:	-20°C/1 year