

Antibodies Product Datasheet

Product Name:	ALS Antibody
Catalog Number:	C20507
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IGFALS. AA range:303-352
Concentration:	1 mg/ml
Gene Name:	IGFALS
Protein Name:	Insulin-like growth factor-binding protein complex acid labile subunit
Other Name:	IGFALS; ALS; Insulin-like growth factor-binding protein complex acid labile subunit; ALS
Human Swiss Prot No.:	P35858
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ALS Polyclonal Antibody detects endogenous levels of ALS protein.
Storage and Stability:	-20°C/1 year