

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | Trk B Antibody |
| Catalog Number: | B0035 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, IHC-p, IF/ICC, ELISA |
| Dilutions: | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications. |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | The antiserum was produced against synthesized peptide derived from human Trk B. AA range:481-530 |
| Concentration: | 1 mg/ml |
| Gene Name: | NTRK2 |
| Protein Name: | BDNF/NT-3 growth factors receptor |
| Other Name: | NTRK2; TRKB; BDNF/NT-3 growth factors receptor; GP145-TrkB; Trk-B; Neurotrophic tyrosine kinase receptor type 2; TrkB tyrosine kinase; Tropomyosin-related kinase B |
| Human Swiss Prot No.: | Q16620 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | Trk B Polyclonal Antibody detects endogenous levels of Trk B protein. |
| Storage and Stability: | -20°C/1 year |