

Antibodies Product Datasheet

Product Name:	INO80B Antibody
Catalog Number:	C49958
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human INO80B
Concentration:	1 mg/ml
Gene Name:	INO80B HMGA1L4 PAPA1 ZNHIT4
Protein Name:	INO80 complex subunit B (High mobility group AT-hook 1-like 4) (IES2 homolog) (hles2) (PAP-1-associated protein 1) (PAPA-1) (Zinc finger HIT domain-containing protein 4)
Other Name:	
Human Swiss Prot No.:	Q9C086
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of INO80B at Human, Mouse
Storage and Stability:	Store at -20° C/1 year