

Antibodies Product Datasheet

Product Name:	DD3 Antibody
Catalog Number:	C14395
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AKR1C3. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	AKR1C3
Protein Name:	Aldo-keto reductase family 1 member C3
Other Name:	AKR1C3; DDH1; HSD17B5; KIAA0119; PGFS; Aldo-keto reductase family 1 member C3; 17-beta-hydroxysteroid dehydrogenase type 5; 17-beta-HSD 5; 3-alpha-HSD type II; brain; 3-alpha-hydroxysteroid dehydrogenase type 2; 3-alpha-HSD type 2; Chlordec
Human Swiss Prot No.:	P42330
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DD3 Polyclonal Antibody detects endogenous levels of DD3 protein.
Storage and Stability:	-20°C/1 year