
Antibodies Product Datasheet

Product Name:	GLU2B Antibody
Catalog Number:	R12-2817
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	GLU2B Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GLU2B.
Concentration:	1 mg/ml
Gene Name:	PRKCSH
Other Name:	80K-H; AGE-R2; AGE-binding receptor 2; G19P1; GLU2B; PCLD; PKCSH; glucosidase II beta subunit; glucosidase II, beta subunit; hepatocystin; protein kinase C substrate; 60.1 kDa protein, heavy chain; P14314
Human Swiss Prot No.:	P14314
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	GLU2B Antibody antibody detects endogenous levels of GLU2B.
Storage and Stability:	Stable at -20°C for at least 1 year.