

Antibodies Product Datasheet

Product Name:	CaMK2 (Phospho-Thr286) Antibody
Catalog Number:	A7033sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA(peptide)1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized pospho peptide derived from human CaMK2 (Phospho-Thr286)
Concentration:	1 mg/ml
Gene Name:	CAMK2D CAMKD
Protein Name:	CaMK2 (Phospho-Thr286)
Other Name:	Calcium/calmodulin-dependent protein kinase type II subunit delta (CaM kinase II subunit delta) (CaMK-II subunit delta) (EC 2.7.11.17)
Human Swiss Prot No.:	Q9UQM7/Q13557
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous pospho levels of human CaMK2 (Phospho-Thr286)
Storage and Stability:	-20°C/1 year