

Antibodies Product Datasheet

Product Name:	Pitx1 Antibody
Catalog Number:	C10791
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PITX1. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	PITX1
Protein Name:	Pituitary homeobox 1
Other Name:	PITX1; BFT; PTX1; Pituitary homeobox 1; Hindlimb-expressed homeobox protein backfoot; Homeobox protein PITX1; Paired-like homeodomain transcription factor 1
Human Swiss Prot No.:	P78337
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Pitx1 Polyclonal Antibody detects endogenous levels of Pitx1 protein.
Storage and Stability:	-20°C/1 year