

Antibodies Product Datasheet

Product Name:	ZIP8 Antibody
Catalog Number:	C48038
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ZIP8
Concentration:	1 mg/ml
Gene Name:	SLC39A8 BIGM103 ZIP8 PP3105
Protein Name:	Zinc transporter ZIP8 (BCG-induced integral membrane protein in monocyte clone 103 protein) (LIV-1 subfamily of ZIP zinc transporter 6) (LZT-Hs6) (Solute carrier family 39 member 8) (Zrt- and Irt-like protein 8) (ZIP-8)
Other Name:	
Human Swiss Prot No.:	Q9C0K1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ZIP8 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year