
Antibodies Product Datasheet

Product Name:	ADRB2 Antibody
Catalog Number:	R12-2457
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	ADRB2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ADRB2.
Concentration:	1 mg/ml
Gene Name:	ADRB2
Other Name:	Beta-2 adrenergic receptor; Beta-2 adrenoceptor; Beta-2 adrenoreceptor; ADRB2; ADRB2R; B2AR;
Human Swiss Prot No.:	P07550
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ADRB2 Antibody antibody detects endogenous levels of ADRB2.
Storage and Stability:	Stable at -20°C for at least 1 year.