

Antibodies Product Datasheet

Product Name:	EMR4 Antibody
Catalog Number:	G246
Reactivity:	Human, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EMR4P. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	EMR4P
Protein Name:	Putative EGF-like module-containing mucin-like hormone receptor-like 4
Other Name:	EMR4P; EMR4; GPR127; PGR16; Putative EGF-like module-containing mucin-like hormone receptor-like 4; EGF-like module receptor 4; G-protein coupled receptor 127; G-protein coupled receptor PGR16
Human Swiss Prot No.:	Q86SQ3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EMR4 Polyclonal Antibody detects endogenous levels of EMR4 protein.
Storage and Stability:	-20°C/1 year