
Antibodies Product Datasheet

Product Name:	FAS Antibody
Catalog Number:	C30400
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human FAS. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	FAS
Protein Name:	Tumor necrosis factor receptor superfamily member 6
Other Name:	FAS; APT1; FAS1; TNFRSF6; Tumor necrosis factor receptor superfamily member 6; Apo-1 antigen; Apoptosis-mediating surface antigen FAS; FASLG receptor; CD95
Human Swiss Prot No.:	P25445
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FAS Polyclonal Antibody detects endogenous levels of FAS protein.
Storage and Stability:	-20°C/1 year