

Antibodies Product Datasheet

Product Name:	KV9.2 Antibody
Catalog Number:	C17827
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KCNS2. AA range:197-246
Concentration:	1 mg/ml
Gene Name:	KCNS2
Protein Name:	Potassium voltage-gated channel subfamily S member 2
Other Name:	KCNS2; KIAA1144; Potassium voltage-gated channel subfamily S member 2; Delayed-rectifier K(+) channel alpha subunit 2; Voltage-gated potassium channel subunit Kv9.2
Human Swiss Prot No.:	Q9ULS6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KV9.2 Polyclonal Antibody detects endogenous levels of KV9.2 protein.
Storage and Stability:	-20°C/1 year