

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Cytokeratin 10 [ABT-CK10]
<b>Catalog Number:</b>	V0142
<b>Clone ID:</b>	ABT-CK10
<b>Reactivity:</b>	Human
<b>Applications:</b>	IHC
<b>Application Dilutions</b>	IHC 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	P13645
<b>Gene Name:</b>	KRT10
<b>Other Name:</b>	Keratin, type I cytoskeletal 10; Cytokeratin-10; CK-10; K10;KPP
<b>Isotype</b>	Mouse IgG
<b>Localization:</b>	N/A
<b>Positive Tissue Control:</b>	N/A
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	Anti-Cytokeratin 10 [ABT-CK10] detects endogenous levels of total Cytokeratin 10 protein.
<b>Storage and Stability:</b>	-20°C/1 year