

Antibodies Product Datasheet

Product Name:	TNF- α Antibody
Catalog Number:	C30026
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human TNF. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	TNF
Protein Name:	Tumor necrosis factor
Other Name:	TNF; TNFA; TNFSF2; Tumor necrosis factor; Cachectin; TNF-alpha; Tumor necrosis factor ligand superfamily member 2; TNF-a
Human Swiss Prot No.:	P01375
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TNF- α Polyclonal Antibody detects endogenous levels of TNF- α protein.
Storage and Stability:	-20°C/1 year