

Antibodies Product Datasheet

Product Name:	Elf-4 Antibody
Catalog Number:	C20675
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ELF4. AA range:178-227
Concentration:	1 mg/ml
Gene Name:	ELF4
Protein Name:	ETS-related transcription factor Elf-4
Other Name:	ELF4; ELFR; MEF; ETS-related transcription factor Elf-4; E74-like factor 4; Myeloid Elf-1-like factor
Human Swiss Prot No.:	Q99607
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Elf-4 Polyclonal Antibody detects endogenous levels of Elf-4 protein.
Storage and Stability:	-20°C/1 year