

Antibodies Product Datasheet

Product Name:	MDA5 Phospho-Ser828 Antibody
Catalog Number:	A9842
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MDA5 ser828
Concentration:	1 mg/ml
Gene Name:	IFIH1 MDA5 RH116
Protein Name:	MDA5 ser828
Other Name:	Interferon-induced helicase C domain-containing protein 1 (EC 3.6.4.13) (Clinically amyopathic dermatomyositis autoantigen 140 kDa) (CADM-140 autoantigen) (Helicase with 2 CARD domains) (Helicard) (Interferon-induced with helicase C domain protein 1) (Mel
Human Swiss Prot No.:	Q9BYX4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.224% sodium azide.
Specificity:	This antibody detects endogenous levels of MDA5 ser828 at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year