
Antibodies Product Datasheet

Product Name:	IGFBP-3 Antibody
Catalog Number:	R12-2914
Reactivity:	Human, Mouse
Applications:	WB, IHC, ELISA
Dilutions:	IGFBP-3 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IGFBP-3.
Concentration:	1 mg/ml
Gene Name:	IGFBP3
Other Name:	IBP-3; IBP3; IGF-binding protein 3; IGFBP3; Insulin-like growth factor binding protein 3 precursor
Human Swiss Prot No.:	P17936
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	IGFBP-3 Antibody antibody detects endogenous levels of IGFBP-3.
Storage and Stability:	Stable at -20°C for at least 1 year.