

## Antibodies Product Datasheet

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | PRMT1 mouse mAb   |
| <b>Catalog Number:</b>        | N0884   |
| <b>Reactivity:</b>            | Human   |
| <b>Applications:</b>          | WB  |
| <b>Dilutions:</b>             | wb 1:1000   |
| <b>Modification:</b>          | Unmodified/ Total   |
| <b>Source:</b>                | Monoclonal Mouse  |
| <b>Purification:</b>          | The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.  |
| <b>Immunogen:</b>             | Purified recombinant human PRMT1 protein fragments expressed in E.coli.   |
| <b>Concentration:</b>         | 1 mg/ml   |
| <b>Gene Name:</b>             | prmt1   |
| <b>Protein Name:</b>          |   |
| <b>Other Name:</b>            | ANM 1;ANM1;ANM1_HUMAN;complete cds;HCP 1;HCP1;Heterogeneous nuclear ribonucleoprotein methyltransferase 1 like 2;Heterogeneous nuclear ribonucleoproteins methyltransferase like 2;Heterogeneous nuclear ribonucleoproteins methyltransferase like2;Histone-arg |
| <b>Human Swiss Prot No.:</b>  | Q99873  |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Specificity:</b>           | This antibody detects endogenous levels of PRMT1 and does not cross-react with related proteins.  |
| <b>Storage and Stability:</b> | -20°C/1 year  |