

Antibodies Product Datasheet

Product Name:	HJURP Antibody
Catalog Number:	C44097
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 80-160
Concentration:	1 mg/ml
Gene Name:	HJURP FAKTS FLEG1 URLC9
Protein Name:	Holliday junction recognition protein (14-3-3-associated AKT substrate) (Fetal liver-expressing gene 1 protein) (Up-regulated in lung cancer 9)
Other Name:	
Human Swiss Prot No.:	Q8NCD3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	HJURP Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year