
Antibodies Product Datasheet

Product Name:	ER α (Acetyl-Lys268) Antibody
Catalog Number:	D12238
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ER α (Acetyl Lys268)
Concentration:	1 mg/ml
Gene Name:	ESR1 ESR NR3A1
Protein Name:	ER α (Acetyl Lys268)
Other Name:	Estrogen receptor (ER) (ER-alpha) (Estradiol receptor) (Nuclear receptor subfamily 3 group A member 1)
Human Swiss Prot No.:	P03372
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human ER α (Acetyl Lys268)
Storage and Stability:	-20°C/1 year