

Antibodies Product Datasheet

Product Name:	MLL Monoclonal Antibody
Catalog Number:	N0443
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of MLL (aa3751-3968) expressed in E. Coli.
Concentration:	
Gene Name:	MLL
Protein Name:	Histone-lysine N-methyltransferase MLL
Other Name:	MLL; ALL1; CXXC7; HRX; HTRX; KMT2A; MLL1; TRX1; Histone-lysine N-methyltransferase MLL; ALL-1; CXXC-type zinc finger protein 7; Lysine N-methyltransferase 2A; KMT2A; Trithorax-like protein; Zinc finger protein HRX
Human Swiss Prot No.:	Q03164
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	MLL Monoclonal Antibody detects endogenous levels of MLL protein.
Storage and Stability:	-20°C/1 year