

Antibodies Product Datasheet

Product Name:	BCKDHB Antibody
Catalog Number:	C49041
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BCKDHB
Concentration:	1 mg/ml
Gene Name:	BCKDHB
Protein Name:	2-oxoisovalerate dehydrogenase subunit beta, mitochondrial (EC 1.2.4.4) (Branched-chain alpha-keto acid dehydrogenase E1 component beta chain) (BCKDE1B) (BCKDH E1-beta)
Other Name:	
Human Swiss Prot No.:	P21953
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BCKDHB at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year