

Antibodies Product Datasheet

Product Name:	TFDP1 Antibody
Catalog Number:	C0175
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DP-1. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	TFDP1
Protein Name:	Transcription factor Dp-1
Other Name:	TFDP1; DP1; Transcription factor Dp-1; DRTF1-polypeptide 1; DRTF1; E2F dimerization partner 1
Human Swiss Prot No.:	Q14186
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TFDP1 Polyclonal Antibody detects endogenous levels of TFDP1 protein.
Storage and Stability:	-20°C/1 year