

Antibodies Product Datasheet

Product Name:	PFK-2 car (Phospho Ser466) Antibody
Catalog Number:	C51147
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human PFKFB2 (Phospho Ser466)
Concentration:	1 mg/ml
Gene Name:	PFKFB2
Protein Name:	6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 2 (6PF-2-K/Fru-2,6-P2ase 2) (PFK/FBPase 2) (6PF-2-K/Fru-2,6-P2ase heart-type isozyme) [Includes: 6-phosphofructo-2-kinase (EC 2.7.1.105)Fructose-2,6-bisphosphatase (EC 3.1.3.46)]
Other Name:	6-phosphofructo-2-kinase/fructose-2,6-bisphosphatase 2 (6PF-2-K/Fru-2,6-P2ase 2) (PFK/FBPase 2) (6PF-2-K/Fru-2,6-P2ase heart-type isozyme) [Includes: 6-phosphofructo-2-kinase (EC 2.7.1.105); Fructose-2,6-bisphosphatase (EC 3.1.3.46)]
Human Swiss Prot No.:	O60825
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PFKFB2 (Phospho Ser466) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year