

---

## Antibodies Product Datasheet

<b>Product Name:</b>	ILK-1 (Phospho-Ser343) Antibody
<b>Catalog Number:</b>	A1069
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human ILK-1 (Phospho-Ser343)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ILK ILK1 ILK2
<b>Protein Name:</b>	ILK-1 (Phospho-Ser343)
<b>Other Name:</b>	Integrin-linked protein kinase (EC 2.7.11.1) (59 kDa serine/threonine-protein kinase) (ILK-1) (ILK-2) (p59ILK)
<b>Human Swiss Prot No.:</b>	Q13418
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.129% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of ILK-1 (Phospho-Ser343) at Human, Mouse,Rat
<b>Storage and Stability:</b>	-20°C/1 year