

Antibodies Product Datasheet

Product Name:	HDAC2 (Phospho-Ser394) Antibody
Catalog Number:	A7099
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HDAC2 around the phosphorylation site of Ser394. AA range:360-409
Concentration:	1 mg/ml
Gene Name:	HDAC2
Protein Name:	Histone deacetylase 2
Other Name:	HDAC2; Histone deacetylase 2; HD2
Human Swiss Prot No.:	Q92769
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-HDAC2 (S394) Polyclonal Antibody detects endogenous levels of HDAC2 protein only when phosphorylated at S394.
Storage and Stability:	-20°C/1 year