

Antibodies Product Datasheet

Product Name:	PR Antibody
Catalog Number:	B0558
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Progesterone Receptor. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	PGR
Protein Name:	Progesterone receptor
Other Name:	PGR; NR3C3; Progesterone receptor; PR; Nuclear receptor subfamily 3 group C member 3
Human Swiss Prot No.:	P06401
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PR Polyclonal Antibody detects endogenous levels of PR protein.
Storage and Stability:	-20°C/1 year