

## Antibodies Product Datasheet

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | MMP19 (Cleaved-Tyr98) Antibody  |
| <b>Catalog Number:</b>        | L0466   |
| <b>Reactivity:</b>            | Human, Mouse  |
| <b>Applications:</b>          | WB, ELISA   |
| <b>Dilutions:</b>             | WB 1:1000-2000 ELISA 1:5000-20000   |
| <b>Modification:</b>          | Cleaved Specific  |
| <b>Source:</b>                | Polyclonal Rabbit   |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                             |
| <b>Immunogen:</b>             | Synthesized peptide derived from human MMP19 (Cleaved-Tyr98)  |
| <b>Concentration:</b>         | 1 mg/ml   |
| <b>Gene Name:</b>             | MMP19 MMP18 RASI  |
| <b>Protein Name:</b>          | MMP19 (Cleaved-Tyr98)   |
| <b>Other Name:</b>            | Matrix metalloproteinase-19 (MMP-19) (EC 3.4.24.-) (Matrix metalloproteinase RASI) (Matrix metalloproteinase-18) (MMP-18)             |
| <b>Human Swiss Prot No.:</b>  | Q99542  |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Specificity:</b>           | This antibody detects endogenous levels of Human, Mouse MMP19 (Cleaved-Tyr98, protein was cleaved amino acid sequence between 97-98 ) |
| <b>Storage and Stability:</b> | -20°C/1 year  |