

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Human Cytokeratin 8 (ABT536) Mouse Monoclonal Antibody
<b>Catalog Number:</b>	V0492
<b>Clone ID:</b>	ABT536
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC, WB
<b>Application Dilutions</b>	IHC 1:100-200, WB 1:500-1000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	KRT8
<b>Gene Name:</b>	P05787
<b>Other Name:</b>	CARD2, CK 8, CK-8, CK8, CYK8, CYKER, Cytokeratin endo A, Cytokeratin-8, DreK8, EndoA, K2C8, K2C8_HUMAN, K8, Keratin 8, Keratin type II cytoskeletal 8, Keratin, type II cytoskeletal 8, Keratin-8, KO, Krt 2.8, KRT8, MGC118110, MGC174782, MGC53564, MGC85764, sb:cb186, Type-II keratin Kb8
<b>Isotype</b>	IgG2a, kappa
<b>Localization:</b>	Liver/ Tonsil
<b>Positive Tissue Control:</b>	Cytoplasmic, Membranous
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	The antibody can specifically recognize human CK8 protein, and shows no cross reaction with CK1, 4, 5, 6, 7, 10, 15, 16, 17, 18, 20.
<b>Storage and Stability:</b>	-20°C/1 year