

Antibodies Product Datasheet

Product Name:	Amylase Antibody
Catalog Number:	C30063
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human AMY1/2. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	AMY1A/AMY1B/AMY1C/AMY2A/AMY2B
Protein Name:	Alpha-amylase 1/Pancreatic alpha-amylase/Alpha-amylase 2B
Other Name:	AMY1A; AMY1; AMY1B; AMY1; AMY1C; AMY1; Alpha-amylase 1; 1,4-alpha-D-glucan glucanohydrolase 1; Salivary alpha-amylase; AMY2A; Pancreatic alpha-amylase; PA;1,4-alpha-D-glucan glucanohydrolase; AMY2B; Alpha-amylase 2B; 1,4-alpha-D-glucan glucanohydrolase 2B
Human Swiss Prot No.:	P04745
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Amylase Polyclonal Antibody detects endogenous levels of Amylase protein.
Storage and Stability:	-20°C/1 year