

Antibodies Product Datasheet

Product Name:	PPIL2 Antibody
Catalog Number:	C49286
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human PPIL2
Concentration:	1 mg/ml
Gene Name:	PPIL2
Protein Name:	Peptidyl-prolyl cis-trans isomerase-like 2 (PPIase) (EC 5.2.1.8) (Cyclophilin-60) (Cyclophilin-like protein Cyp-60) (Rotamase PPIL2)
Other Name:	
Human Swiss Prot No.:	Q13356
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of PPIL2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year