

Antibodies Product Datasheet

Product Name:	SMAD5(C-term) mouse mAb
Catalog Number:	N0824
Reactivity:	Human
Applications:	WB, ICC, FC
Dilutions:	wb 1:1000 icc 1:75 fcm 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human SMAD5 (C-terminus) protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	smad5
Protein Name:	
Other Name:	DKFZp781C1895;DKFZp781O1323;Dwfc;hSmad 5;hSmad5;JV5 1;JV5-1;MAD homolog 5;MAD mothers against decapentaplegic homolog 5;MAD, mothers against decapentaplegic homolog 5;MADH 5;MADH5;Mothers against decapentaplegic homolog 5;Mothers against DPP homolog 5;Mus
Human Swiss Prot No.:	Q99717
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SMAD5 (C-terminus) and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year