

Antibodies Product Datasheet

| | |
|---------------------------------|--|
| Product Name: | Anti-Human Cytokeratin 20 (ABT044) Mouse Monoclonal Antibody |
| Catalog Number: | V0305 |
| Clone ID: | ABT044 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | IHC, WB |
| Application Dilutions | IHC 1:100-200, WB 1:500-1000 |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| UniProt: | KRT20 |
| Gene Name: | P35900 |
| Other Name: | CD20, CK 20, CK-20, CK20, Cytokeratin-20, Cytokeratin20, K1C20_HUMAN, K20, KA20, Keratin 20, keratin 20, type I, keratin 21, rat, homolog of, Keratin, Keratin type I cytoskeletal 20, Keratin-20, Keratin20, KRT 20, KRT 21, KRT20, KRT21, MGC35423, OTTHUMP00000164518, Protein IT, type I cytoskeletal 20 |
| Isotype | IgG2b, Kappa |
| Localization: | Colon/ Appendix |
| Positive Tissue Control: | Cytoplasmic, Membranous |
| Immunogen | Synthetic peptide |
| Purification | Immunogen affinity purified |
| Formulation: | PBS, pH7.4, 50% glycerol, 0.03%Proclin300 |
| Specificity: | The antibody can specifically recognize human CK20 protein, and shows no cross reaction with CK1, 5, 6, 7, 8, 13, 14, 15, 17, 18. |
| Storage and Stability: | -20°C/1 year |