

Antibodies Product Datasheet

Product Name:	SLC25A20 Antibody
Catalog Number:	C48504
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SLC25A20
Concentration:	1 mg/ml
Gene Name:	SLC25A20 CAC CACT
Protein Name:	Mitochondrial carnitine/acylcarnitine carrier protein (Carnitine/acylcarnitine translocase) (CAC) (Solute carrier family 25 member 20)
Other Name:	
Human Swiss Prot No.:	O43772
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SLC25A20 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year